

## RACKS



One of the biggest advantages of the Lancer range is basket interchangeability.

From the 810 LX to the 1300 LX as well as the 1800 LXA, all use the same range of baskets. You gain flexibility while saving money.

From basic baskets to customized ones – Lancer has the right solution to meet your washing needs.

## CONSUMABLES



Lancer LLL detergent and Lancer LLL neutralizing agent are tested and validated – the perfect match for our washers.

## OPTIONS

Our Labexia range offers a wide range of options, from printers and steam condensers to validation of the washing process.

# LANCER

## GETINGE GROUP

Getinge Lancer  
30 Boulevard de l'Industrie  
F-31170 Tournefeuille, France  
Tél: +33 5 61 15 11 11  
Fax: +33 5 61 15 16 16  
lancersite@getinge.com www.lancer.fr

LAB GLASSWARE  
WASHER-DRYERS  
YOUR PARTNER IN  
LABORATORY WASHING

# LANCER

2231rev2 ENG 10.04 - XT - We reserve the right to make technical and construction changes.



## UNDERBENCH WASHERS

When you need to optimize your lab space with no compromise on your professional washing.

### Lancer 810 LX

Start-up kit for professional washing results.

- 2 washing levels
- High loading flexibility
- Economical
- Full professional design
- IQ/OQ
- No drying
- Total chamber volume 138 L



### Lancer 815 LX

Hot forced-air drying system outside the glassware.

- 2 washing levels
- High loading flexibility
- Professional design
- IQ/OQ
- Total chamber volume 138 L



### Lancer 820 LX

Drying inside and outside the glassware, on board storage chemicals.

- 2 washing levels
- High loading flexibility
- HEPA jet drying system
- IQ/OQ
- Total chamber volume 138 L



Our washers are designed and manufactured in an environmentally sound facility certified since 2003 in accordance with ISO140001. Thanks to constant improvement of our products, we have achieved one of the lowest water and energy consumption rates in the industry.

## VERTICAL WASHERS

Large washing capacity with a small footprint.

### Lancer 910 LX

The first step into higher capacity washing.

- 2 washing levels
- Small footprint
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 148 L



### Lancer 1300 LX

The cost-effective solution for high-capacity loads.

- 3 washing levels
- Small footprint
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 175 L



### Lancer 1400 LX

Larger and wider chamber to handle mixed loads.

- 3 washing levels
- Increased capacity
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 296 L



## PREMIUM MODELS

Highest capacity models with increased specifications for more demanding applications.

### Lancer 1400 LXP

Optimal solution to validate your washing cycle.

- 3 washing levels
- Increased capacity
- QA pharmaceutical
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 296 L



### Lancer 1600 LXP

Wash and dry on 4 levels at the same time.

- 4 washing levels
- Increased capacity
- QA pharmaceutical
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 369 L



### Lancer 1800 LXA

Up to 6 separate baskets in 8 different positions.

- 6 washing levels
- Automatic glass door
- QA pharmaceutical
- IQ/OQ
- HEPA jet drying system
- Total chamber volume 417 L

